

A C C U R A

Medizintechnik GmbH

Max-Planck-Str. 33 61184 Karben Germany
Tel +49-6039-9201-0 Fax +49-6039-9201-22
E-Mail info@accura.info
www.accuramed.de



Guidewires

We supply a large range of different wires for medical applications. Our wires show high grade quality having proper **flexibility**, high **tension** and best **torque**.

All wires can be supplied sterile packed in dispenser, OEM or bulk.

Please order according to following specifications:

- Straight tip (S) or
- J-Tip (J) with diameter (mm)
- Fixed core (F) or movable core (M)
- Coated (C) or uncoated (U)
- Length (cm)
- Diameter (inch)



Führungsdrähte

Wir liefern ein sehr umfangreiches Programm medizinischer Führungsdrähte. Sie zeichnen sich insbesondere durch hohe **Zugfestigkeit**, angepaßte **Flexibilität** und gute **Drehstabilität** aus.

Alle Drähte können nach Kundenwunsch steril verpackt im Dispenser oder unsteril geliefert werden.

Bitte bestellen Sie entsprechend nachfolgenden Spezifikationen:

Bestellinformationen:

- Gerade (S) oder
- J-Spitze (J) mit Durchmesser in mm
- Feste Seele (F) oder bewegliche Seele (M)
- Beschichtet (C) oder unbeschichtet (U)
- Länge
- Durchmesser (Inch)