



ACCURA

PERSOON

Peripheral Balloon Dilation Catheter

*Your personal balloon for **quick peripheral revascularisation***

Accura products for a blooming health



Developed under German Technology
Manufactured in Germany



PERSOON - Peripheral Balloon Dilation Catheter

Your personal balloon for quick peripheral revascularisation

- **5-fold balloon** provides a **small profile**
- **Minimal winging**
- **Good radiopaque visibility**
- Available in **4F and 5F**
- **Semi compliant** balloon characteristic
- **Low profile**
- **Flexible** catheter
- **Good pushability**

| | Available Balloon Diameters (mm) | Available Balloon Lengths (cm) | Available Catheter Lengths (cm) | Shaft Diameter (mm) | | |
|----------------------------------|----------------------------------|--------------------------------|---------------------------------|---------------------|---------------------------|--|
| | 3.0 | 1.5 | 40 | 1.35 | | |
| | 4.0 | 2.0 | 50 | 1.67 | | |
| | 5.0 | 3.0 | 75 | | | |
| | 6.0 | 4.0 | 120 | | | |
| | 7.0 | 5.0 | 135 | | | |
| | 8.0 | 6.0 | 150 | | | |
| | 9.0 | 8.0 | | | | |
| | 10.0 | 10.0 | | | | |
| | 12.0 | 12.0 | | | | |
| | 14.0 | | | | | |
| Examples for Reference Code: | | | | Reference Code | | |
| The item number is combined of : | | | | | | |
| The letters "PB" in the front | Balloon Diameter | Balloon Length | Catheter Length | Shaft Diameter | The letter "G" in the end | |
| PB | 7.0 | 4.0 | 120 | 1.67 | PB70401216G | |
| PB | 3.0 | 6.0 | 75 | 1.35 | PB30607513G | |