



ACCURA

# INTRODUCER SET

The smooth Access

*Let us introduce this set to you*

Accura products for a blooming health





ACCURA

## ACCURA INTRODUCER-SET

The Set contains all (or part) of the following components:

- **Dilator** - 19 mm HDPE Dilator with special shaped and tapered tip
- **Sheath** - Sheath, 11 cm / 23 cm length, smoothly adapted to the Dilator and made of FEP
- **Haemostatic Valve** - Haemostatic Valve with side-link, automatic sealing (no screwing needed to fix the catheter)
- **Guide-Wire** - J-Guide-Wire adapted to the Dilator with straightener
- **Seldinger Needle** 18G, 8 cm length
- **Scalpel**

All sets are sterile packed in a transparent tray and are designed for single use.

Available Sizes		
	Sheath length	Item No.
4 F	9cm	3004009
5 F	11cm	3005011
6 F	11cm	3006011
7 F	11cm	3007011
8 F	11cm	3008011
9 F	11cm	3009011
10F	11cm	3010011

Available Sizes (23 cm)		
	Sheath length	Item No.
6 F	23cm	3006023
7 F	23cm	3007023
8 F	23cm	3008023
9 F	23cm	3009023
4 F (Transradial)	9cm	3004009R
5 F (Transradial)	11cm	3005011R
6 F (Transradial)	11cm	3006011R
7 F (Transradial)	11cm	3007011R

Custom made sets with other sizes (lengths) or additional items or fewer components are available upon request as well as Split Introducers. We also offer Dilators 4 F up to 20 F.